

Physics


Chemistry · Biology

Technology



Leybold Didactic GmbH  
Lehr- und Didaktiksysteme

05/02-W97-Hke



## Instruction sheet 524 060

Force sensor S,  $\pm 1\text{N}$  (524 060)

- 1 Holder
- 2 4 mm bore for experiment material
- 3 Experiment supporting points, electrically insulated
- 4 Stand rod
- 5 Plug for Sensor-CASSY


## 1 Description

The force sensor S,  $\pm 1\text{N}$  is intended for use in conjunction with the Sensor-CASSY, to which it can be directly connected. The sensor enables forces in the range  $\pm 1\text{ N}$  to be measured; in all measuring ranges, a tare of up to  $\pm 2.5\text{ N}$  can be compensated for the zero adjustment (range of compensation). As there is no mechanical damping, oscillations can also be investigated.

Experiment descriptions for using the force sensor S are found in the CASSY Lab help.

## 2 Principle of operation

The force sensor consists of a very precise double flexion element with four strain gauges (**a**) in bridge configuration. It can be used in two directions. The flexion element is located in a robust housing, which also limits the maximum mechanical deflection thus providing protection against damage. The required electronics is placed in the plug.



## 3 Scope of supply

- 1 force sensor S,  $\pm 1\text{N}$
- 1 fixing bracket on 4-mm pin

## 4 Technical data


Load	max. $\pm 2.5\text{ N}$
Measuring ranges	$\pm 10\text{ mN}$ , $\pm 100\text{ mN}$ , $\pm 1\text{ N}$
Range of compensation	$\pm 2.5\text{ N}$ in all measuring ranges
Mechanical deflection:	$\pm 0.5\text{ mm/N}$
Measuring error:	$< 1\%$
Resolution:	$< 0.01\text{ mN}$

The useful resolution depends on environmental vibration, draught, temperature variation, etc. Therefore computers with a fan should not be put on the experiment table.

Connections:

SUB D 15 plug with integrated amplifier for connecting to the Sensor-CASSY

## 5 Operation



The force sensor is fixed via a short stand rod and used in a way that forces can be easily measured in two directions (in contrast to electronic balances).

At the free end at the front of the force sensor, there is a plastic part for attaching experiment material such as conductor loops or test balls; arrows indicate the direction in which forces can be measured.

additionally required:

1 Sensor-CASSY

524 010

1 CASSY Lab

524 200

1 PC with Windows

or

1 CASSY-Display

524 020