
Chemie · Biologie

LEYBOLD DIDACTIC GMBH

3/94-Sf-

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The multi-purpose microphone is a battery-operated microphone for sound waves within the frequency range of 30 Hz to 20 kHz and, with reduced sensitivity, also for the ultrasonic range up to 40 kHz. An integrated preamplifier enables direct connection of measuring instruments (e. g. oscilloscope, stop-clocks and counters). An additional level output is provided for direct connection to a chart recorder.

In view of the built-in automatic battery cut-out and the very low power consumption of the preamplifier the battery life is mainly dependent on the shelf life of the battery itself (1 to 2 years).

Examples of experiments:

- Music and speech recordings
- Reflection and interference of sound waves
- Operation of electric timers and frequency counters for the measurement of sound velocity or acoustic frequencies

1 Safety notes

Prevent foreign particles or humidity from entering the tube of the microphone ①.

Max. permissible external voltage applied to connecting

lead: 10 V!

2 Standard equipment, description, technical data (Fig. 1)

- Tube, dia. 8 mm, 25 cm long, mountable into any desired clamp
- ② Mouthpiece with flush microphone capsule
- (3) Handle with preamplifier electronics and battery Dimensions: 19 cm x 3.5 cm x 3.5 cm
- (4) ON pushbutton:
 - By briefly pressing this button the preamplifier is switched on. This is also possible during operation.
 - The device switches off automatically approx. 30 minutes after last switch-on.
- ⑤ Gain adjustment, continuously variable up to a factor of 16

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Position =: ("level"), if a chart recorder is connected.

Position Π : ("trigger"), if any of the following instruments is connected: electronic timers or frequency meters, e. g. LH digital counter (575 40), electronic stop-clock D (313 01), electronic stop-clock P (313 031), counter P (575 45), or digital input of an interface (e.g. from 524 007).

- Output cable, 2 m long, with one black ("earth") and one yellow ("signal") 4-mm plug.
- 8 Battery compartment (for 9-V block battery), with sliding cover
- Stand rod with thread, 22 cm x 1 cm dia., to be screwed into the matching socket in the bottom of the handle.

Not shown: Battery 9 V (according to IEC 6F 22)

Technical data:

Transducer: electret capsule, 5 mm dia., at the tip of the tube.

Data of microphone including preamplifier:

Frequency range: 30 Hz to 20 kHz (40 kHz)

Characteristic: omnidirectional Sensitivity, related to signal output:

approx. 18 to 300 mV/µbar at 1 kHz, depending on gain

Sound-level limit: 100 dB

Signal output:

Operating limit: 4 V_{pp} for connection to instruments with $R_i \ge 20 \text{ k}\Omega$ Output impedance: < 1 k Ω (f > 60 Hz)

Level output:

Output impedance: $10 \text{ k}\Omega$ Operating limit: > 1 VNonlinearity: approx. 50 mV

Trigger output:

Open-collector output Bias voltage: 10 V max. Current: 20 mA max.

Power supply:

9 V block battery (685 45) according to IEC 6F 22

Current drain: 1.3 mA

Dimensions without stand rod: 45 cm x 3.5 cm x 3.5 cm

Weight: 0.45 kg

3 Operation

3.1 Connection to AC voltmeters, oscilloscope or amplifier

Set function switch 6 to \sim ("signal") and select a suitable AC voltage range on measuring instrument. Input impedance of measuring instrument used should exceed > 20 k Ω , otherwise reduction of operating limit (signal limiting).

3.2 Connection to chart recorder

Set function switch 6 to = ("level"). Input impedance of chart recorder \geq 100 k Ω .

3.3 Connection to counter or stop-clock

Set function switch 6 to $\Pi("trigger")$. Always connect black plug of connecting lead to the instrument socket marked with the earth symbol \bot .

3.4 Replacing the battery

Open battery compartment ® by sliding the cover to the side, remove battery, insert new battery (685 45) and close battery compartment.

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Supplementary Instructions for

1 Safety notes

Strong electromagnetic fields can interfere with the sensitive electronics of the multi-purpose microphone. This can result in a temporary impairment of function (e.g. distotion or superpositioning of the desired signal).

Measures:

- Do not operate any RF-transmitters (e.g. cellular phones) in the experiment room which are not a part of the experiment setup.
- Keep all connecting leads which could act as antennas as short as possible. In particular, do not extend the connection lead of the multi-purpose microphone.

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